

Required Materials

| • | Trilogy Laboratory Fluorometer | P/N: 7200-000/002 |
|---|---|-------------------|
| • | Chlorophyll a Acidification Module | P/N: 7200-040(-W) |
| • | 12x75mm round bottom glass test tubes | P/N: 10-029A |
| • | 12mm round vial adaptor (included with Trilogy accessories) | P/N: 016-0810 |
| • | Glass Cuvette – only for use with P/N 7200-040-W | P/N 7000-955 |
| • | Liquid Chlorophyll a Standard – available from Turner Designs | P/N: 10-850 |
| | | |

Optional Materials

Solid Secondary Standard (refer to User Manual for instruction) P/N: 8000-952

Procedure

- 1. Lift the Trilogy Laboratory Fluorometer's lid.
- 2. Snap the Chlorophyll a Acidification Module into the Trilogy Laboratory Fluorometer.
- 3. Insert the 12mm round vial adapter into the module if using glass test tubes.
- 4. Close the lid and turn on the fluorometer.
- 5. Press the "Chl-A" button on the touch screen to choose the acidification module.
- 6. Press "OK" after confirming the correct module is snapped in to go to the home screen.
- 7. From the home screen choose the "Calibrate" button.
- 8. Select "Run New Calibration" and choose "µg/L".
- 9. Insert a 90% acetone blank solution and press "OK". Wait for measuring to complete.
- 10. Remove the blank and insert a standard of known concentration.

Note: If you will enter more than one standard you must enter them in order of increasing concentration and we recommend the first standard entered to be no less than 1µg/L.

- 11. Enter the concentration using the keypad on the touch screen and press "OK".
- 12. Press "OK" to measure before acidification. Wait for measuring to complete.
- 13. Acidify the standard by adding 0.03 mL of 0.1 N HCl for every mL of standard, wait 90 seconds while gently inverting the vial about 10 times and insert into the instrument.
- 14. Press "OK" to measure after acidification. Wait for measuring to complete.
- 15. Acid Ratio will be displayed on the screen, press "OK" to proceed.

Note: An acceptable acid ratio is >1.7; An optimal acid ratio is 1.9-2.0.

If your acid ratio is <1.7, then you can try one of the following three steps to obtain an acceptable acid ratio:

- a. Remake the standard concentration used as your first standard in the calibration
- b. Use a higher standard concentration as your first standard in the calibration (up to 5 μ g/L)
- c. Obtain a new parent standard to remake all of your standard concentrations in the calibration
- d. Remake the 0.1 N HCL acid
- 16. You may enter up to 5 standards increasing in concentration, or proceed with current calibration.

Note: Acid ratio is only calculated for the first standard entered.

17. Save the calibration.





- 18. Extract all filtered samples as described in EPA Method 445.0, rev. 1.2 (September 1997), section 11.1.
- 19. Insert a sample and press the "Measure Fluorescence" button.
- 20. Enter the volume filtered and extraction volume when prompted.
- 21. Press "OK" to measure the fluorescence before acidification.
- 22. Acidify the sample by adding 0.03 mL of 0.1 N HCl for every mL of sample, wait 90 seconds while gently inverting the vial about 10 times and insert into the instrument.
- 23. Press "OK" to measure the fluorescence after acidification.
- 24. The corrected chlorophyll *a* and pheophytin *a* concentrations in the whole water sample, according to <u>EPA</u> <u>Method 445.0</u> equations in sections 12.2 and 12.3, will be displayed on the home screen.

