

Congratulations on your new *AquaFluor* Handheld Fluorometer/Turbidimeter. We are committed to customer satisfaction. If you need assistance, technical specialists are available to answer your questions at 408-749-0994 or toll-free at 877-316-8049. This Quick Start Guide will help you set up the *AquaFluor* Handheld Fluorometer/Turbidimeter to begin collecting data as quickly as possible.



1 Setup

Before the *AquaFluor* can be used, the supplied batteries must be installed.

1. On the back side of the instrument, remove the battery panel by loosening the center screw and then gently pulling on the screw.
2. Install the 4 - AAA batteries and verify the batteries are positioned correctly.
3. Replace the battery panel and press it down into position.
4. Hold the battery panel down in place and gently tighten the screw until it is snug. Do not over tighten the screw.

Note: If the battery panel is very difficult to press down or remove, you can apply a small amount of silicon-based grease to lubricate the o-ring as needed.

2 Calibration

1. Turn the *AquaFluor* on by pressing <ON/OFF>. Wait 5 seconds for the *AquaFluor* to warm up.
2. Press <A/B> to toggle between different channels and select the channel you want to use. To identify the configuration look at the label on the back side.
3. Assign a Calibration Standard Value. This is the numeric value that you want the standard to read.
4. Press the <STD VAL> button. Use the ↑ and ↓ arrow buttons to set the standard value.
5. To perform a calibration, press the <CAL> button.
6. Press <ENT> to start the calibration.
7. Insert your blank sample and press <ENT>. *AquaFluor* will average the reading for 10 seconds and set the blanking zero point.
8. Insert the standard sample and press <ENT>. The reading is averaged for 10 seconds and the Standard Calibration value is set.
9. Press <ENT> when the calibration is complete to accept the calibration.
10. Press <DIAG> and wait for %FS-BLK and %FS-STD to be displayed. The ratio of the %FS-STD:%FS-BLK should be >5.
11. If the calculate value is not >5, you should use a more concentrated standard to achieve this value.
12. If you have more than one channel you can press <A/B> to toggle to the other channel and repeat steps 3-10.
13. You are now ready to begin analyzing samples.

3 Sample Analysis

1. Insert your sample. The *AquaFluor* is designed with “ambient light rejection”. The black sample compartment cover does NOT need to be closed when reading samples.
2. The orientation and cleanliness of the cuvettes can have an impact on the accuracy of your results. Use clean cuvettes which are dry on the outside. Mark your cuvettes for consistent orientation. Refer to Section 5 in the Users Manual for additional details.
3. Press the <READ> button. The instrument will measure and average the fluorescence signal for 10 seconds.
4. The reading result will be displayed on the top line of the Home screen.
5. The top left corner will then display “WAIT” for 3 seconds. Once “WAIT” disappears, another sample reading can be performed.