



Growth Responses in the Marine Picocyanobacteria *Synechococcus*: A Multi Stressor Approach



Samarpita Basu* and **Katherine R. M. Mackey**
 Department of Earth System Science, University of California Irvine
 *samarpb@uci.edu



Background

- Marine picocyanobacteria are the most widely distributed group of photosynthetic prokaryotes on Earth, and contribute up to 50% of fixed C in the open ocean¹.
- It is currently unknown how global change such as ocean acidification will impact picocyanobacteria on a global scale.
- Carbon dioxide (CO₂) is one of the principal drivers of global change².
- While CO₂ concentration in the atmosphere is estimated to be about 270 ppm before the industrial revolution, it has currently increased to about **415 ppm** and is expected to reach **800–1000 ppm** by the end of this century according to the “business as usual” CO₂ emission scenario³.
- Marine ecosystems are a major sink for atmospheric CO₂, currently accounting for the removal of one-quarter to one-third of anthropogenic CO₂ emissions from the atmosphere⁴.
- The dissolution of anthropogenic CO₂ in the ocean and the subsequent formation of carbonic acid has already resulted in a 30% (0.1 pH unit) increase in seawater [H⁺] and will continue to lower pH by an additional 0.2-0.3 pH units by the end of the century. This decline in ocean pH is referred to as **ocean acidification**⁵.
- At the same time, warming will increase the mean surface temperatures by an average of 3 °C, leading to longer periods of **stratification** with fewer deep mixing events³ (Fig. 1).
- To date, most ocean acidification studies have focused on large-celled, bloom forming phytoplankton such as diatoms and coccolithophores, and there are **very few reports on picocyanobacteria**⁶.




Figure 1. Global change effects on the surface ocean: By the year 2100, pH of the ocean will decline to 7.8 due to increased uptake of atmospheric CO₂. Concomitantly, increased thermal stratification will trap phytoplankton in the surface ocean, resulting in increased light exposure and lower nutrient availability to the cells³.

Methods

- The *Synechococcus* strain **WH 8109** (previously isolated from northwestern Atlantic Ocean slope water) was used in this study.
- Cultures of WH 8109 were maintained in dilute, semi-continuous batch culture in a temperature controlled growth chamber (22 °C) under continuous white light.
- The culture media consisted of 75% coastal seawater and 25% MilliQ water with SN nutrient amendments⁷.
- Cultures were grown at three different temperatures (22 °C, 24 °C and 26 °C) and three CO₂ levels (400 ppm, 600 ppm and 800 ppm) simultaneously (Fig. 2a).
- A constant light intensity of 20 μE m⁻² s⁻¹ was maintained during the experiments.
- Growth was monitored via measurement of **raw fluorescence** on the Turner fluorometer (**Turner Designs Trilogy**).

$$\mu = \frac{\ln(X_1/X_0)}{t_1 - t_0}$$

(Where X₁ and X₀ are the fluorescence values obtained on the days t₁ and t₀ respectively).

- Specific growth rates μ (d⁻¹) were calculated using two time points from within the exponential phase of the growth curve using the following equation:
- The photochemical efficiency of PS II in the dark adapted state (F_v/F_m) and the functional absorption cross section of PSII in the dark adapted state (σ_{PSII}) were monitored on a FIRE fluorometer.
- The daily changes in the medium pH were also monitored.





Figure 2a. Experimental set up consisting of 3 glass tanks with temperature controlled water for culturing *Synechococcus* WH8109 under 3 different temperature and CO₂ levels.

Figure 2b. Raw fluorescence measurement conducted on Turner Designs Trilogy Laboratory Fluorometer.

Results



Specific growth rates of *Synechococcus* WH8109 incubated at 22 °C (Fig. 3a), 24 °C (Fig. 3b) and 26 °C (Fig. 3c) under varying CO₂ concentrations.

F_v/F_m values of *Synechococcus* WH8109 incubated at 22 °C (Fig. 4a), 24 °C (Fig. 5a) and 26 °C (Fig. 6a) under varying CO₂ concentrations, and σ_{PSII} values of *Synechococcus* WH8109 incubated at 22 °C (Fig. 4b), 24 °C (Fig. 5b) and 26 °C (Fig. 6b) under varying CO₂ concentrations.

Conclusions

- The highest specific growth rate of 0.47 d⁻¹ was observed at 22 °C and 800 ppm CO₂ concentration.
- Doubled CO₂ concentration of 800 ppm resulted in higher specific growth rates in *Synechococcus* strain WH 8109 under all the three temperature treatments.
- The values of photochemical efficiency of PS II in the dark adapted state (F_v/F_m) for *Synechococcus* strain WH 8109 remained at or near maximal levels at 26 °C.
- Therefore, increased CO₂ availability combined with higher temperatures may benefit the picocyanobacteria *Synechococcus*, resulting in their dominance in the future ocean.
- However, because individual picocyanobacteria strains may respond differently to future CO₂ and temperature increases, further investigation with additional strains is required prior to generalizing their responses to global change.

Acknowledgements

This work was supported by a Clare Boothe Luce endowment, an Alfred P Sloan Research Fellowship in Ocean Sciences, the Marion Milligan Mason Award for Women in the Chemical Sciences from AAAS and a Simons Early Career Award in Marine Microbial Ecology and Evolution to KRMM. We acknowledge the sponsorship provided by Turner Designs to present in the Ocean Sciences Meeting 2020.



References

- Fu, F. X., Warner, M. E., Zhang, Y., Feng, Y., & Hutchins, D. A. Effects of increased temperature and CO₂ on photosynthesis, growth, and elemental ratios in marine *Synechococcus* (cyanobacteria). *Journal of Phycology* **2007**, 43(3), 485-496.
- Lam, M.K.; Lee, K.T.; Mohamed, A.R. Current status and challenges on microalgae-based carbon capture. *International Journal of Greenhouse Gas Control* **2012**, 10, 456-469.
- Basu, S.; Mackey, K. R. M. Phytoplankton as key mediators of the biological carbon pump: Their responses to a changing climate. *Sustainability* **2018**, 10, 869.
- Li, W.; Gao, K.S.; Beardall, J. Interactive effects of ocean acidification and nitrogen-limitation on the diatom *Phaeodactylum tricornutum*. *PLoS ONE* **2012**, 7, e51590.
- Bhadury, P. Effects of ocean acidification on marine invertebrates-a review. *Indian Journal of Geo-Marine Sciences* **2015**, 44(4), 1-11.
- Mackey, K. R., Morris, J. J., Morel, F. M., & Kranz, S. A. Response of photosynthesis to ocean acidification. *Oceanography* **2015**, 28(2), 74-91.
- Waterbury, J.B., Watson, S.W., Valois, F.W., and Franks, D.G. Biological and ecological characterization of the marine unicellular cyanobacterium *Synechococcus*. *Can Bull Fish Aquat Sci* **1986**, 214: 120.