

10AU Field Fluorometer



10AU Field Fluorometer

The 10AU Field Fluorometer is a rugged, field-portable instrument that can be set up for continuous-flow monitoring or discrete sample analyses. Measurements are extremely stable with stability measured in months or weeks as opposed to days or hours. Watertight filter paddles enable quick change of excitation and emission filters so a variety of compounds can be easily measured on-site.

Available Optical Kits

Optical kits with appropriate filters and lamps are available for a wide range of applications and can be easily installed in the field or lab. Individual filters and lamps can also be purchased for customized applications.

Available Optical Kits

- Ammonium
- Blue Green Algae - Phycocyanin (freshwater)
- Blue Green Algae - Phycoerythrin (marine)
- CDOM/FDOM
- Chlorophyll *a* - Extracted and/or *In Vivo*
- Crude Oil (Long Wavelength)
- Fluorescein Dye Tracing
- Histamine
- Optical Brighteners for Wastewater Monitoring
- Refined Oil (Short Wavelength)
- Rhodamine WT Dye Tracing
- Turbidity

Contact us for Custom Optical Kits

10AU Highlights

- Exceptional Sensitivity
- Automatic or Manual Ranging
- Analog and Digital Output
- Field Installable Optical Kits
- Easy to Read LCD Display

User Configured Options

- Standard or Red PMT
- Continuous Flow Monitoring or Discrete Sample Analysis



Optional Accessories

- Internal Data Logging & Electronic Chart Recording
- Temperature Compensation for Continuous Flow Measurements
- Solid Secondary Standard

Reliable Instruments for an Unreliable World

10AU Field Fluorometer



Ordering Information

BASE INSTRUMENT W/ WATERTIGHT CASE

| | PART NUMBER |
|------------------------------------|-------------|
| 10AU Fluorometer with Red PMT | 1100-100 |
| 10AU Fluorometer with Standard PMT | 1100-101 |

SAMPLE CELL

| | PART NUMBER |
|--|-------------|
| 25 mm One-Piece Continuous Flow System | 1100-020 |
| 25 mm & 13 mm Discrete Sample Test Tube Holder Set | 10-AU-030 |

SOFTWARE

| | PART NUMBER |
|--|-------------|
| Internal Data Logging & Electronic Chart Recording | 10-AU-450 |
| Temperature Compensation Package (On-Line Systems) | 1100-500 |

OPTICAL KITS

| | PART NUMBER |
|---|-------------|
| Ammonium/CDOM/FDOM | 10-303 |
| Blue Green Algae - Phycocyanin (freshwater) | 10-305 |
| Blue Green Algae - Phycoerythrin (marine) | 10-304 |
| Chlorophyll <i>a</i> (<i>In Vivo</i> & Extractive) | 10-037R |
| Chlorophyll <i>a</i> (Extractive) | 10-040R |
| Chlorophyll (<i>In Vivo</i>) | 10-096R |
| Crude Oil (Long Wavelength Oil) | 10-302R |
| Fluorescein Dye | 10-086R |
| Optical Brighteners for Wastewater Monitoring | 10-306 |
| Refined Oil (Short Wavelength Oil) | 10-301R |
| Rhodamine WT Dye | 10-041R |
| Turbidity | 10-307R |

Contact us for Custom Optical Kits

OPTIONAL ACCESSORIES

| | PART NUMBER |
|---|-------------|
| Solid Secondary Standard | 10-AU-904 |
| Shipping Case | 10-AU-060 |
| 12 VDC Power & Signal Cable w/ Clips | 10AU-810 |
| Universal Power & Signal Cable, 115 VAC & 230 VAC | 10AU-800 |

Contact Us

Toll-Free : 1.877.316.8049

Email : sales@turnerdesigns.com

Address:

Phone : 408.749.0994

Web : www.turnerdesigns.com

1995 N 1st Steet

Fax : 408.749.0998

San Jose, CA 95112

10AU Electrical Specifications

| | |
|----------------|---|
| Light Source | Lamp (4 watts, 8000 hours lamp life) |
| Detector | PhotoMultiplier Tube; Standard (300-650 nm), Red (185-870 nm) |
| AC Power | 100-130 V; 200-240 V, 50/60 Hz, 30 watts |
| DC Power | 11-16 V; 2.5 amperes |
| Digital Output | 100% ASCII data output via RS-232 at 4800 or 9600 bps |
| Analog Output | 0.1, 1, 2, or 5 volts |

10AU Detection Limits

| | |
|--------------------------------|---------------------------------------|
| Extracted Chlorophyll <i>a</i> | 0.025 µg/L using EPA Method 445.0 |
| Rhodamine WT Dye | 0.01 ppb RWT in potable water |
| Fluorescein Dye | 0.01 ppb Fluorescein in potable water |

10AU Linear Range

| | |
|--------------------------------|--------------|
| Extracted Chlorophyll <i>a</i> | 0 - 250 µg/L |
| Rhodamine WT Dye | 0 - 250 ppb |
| Fluorescein Dye | 0 - 250 ppb |

10AU Physical Specifications

| | |
|-----------------------|------------------------------|
| Weight | 12.6 kg (28.5 lbs) |
| Height | 24 cm (9.45 in) |
| Width | 55 cm (21.65 in) |
| Depth | 34 cm (13.39 in) |
| Operating Temperature | 0-55° C, 32-131° F (ambient) |