

Certificate of Analysis

Spectrophotometric Chlorophyll Standards

12 July 2021

Lot CAS 167-03

Original Chlorophyll Solution A:

Mixed: 30 March 2021

Spectrophotometric Data: 3 June 2021

$Abs_{664} = 0.87366$ $Abs_{750} = -8.5602E-3$

Chlorophyll concentration determination (Jeffrey Method, 1997)

Chlorophyll-a = 10.6 mg/L

Chlorophyll dilution A is provided with this note. The concentration value is guaranteed through June, 2022, but degradation depends on the handling of the standard. Standards must be kept at temperatures below -20°C and in the dark until use. Spectrophotometers should be standardized in subdued light with this solution warmed to room temperature, about 22°C . The determination of chlorophyll-a in aqueous acetone solvent is a sensitive technique. Please pay special attention to the configuration of the spectrophotometer itself. A bandpass of 2.0 nm or less is necessary to properly use this standard.



John Salinas
Analytical Chemist

12 July 2021
Date

TURNER
 **DESIGNS**

Reliable Instruments for an Unreliable World
1995 N. 1st Street Telephone 408-749-0994
San Jose, CA 95112 Toll Free 877-316-8049
www.turnerdesigns.com Fax 408-749-0998